



DDC Dolphin Ltd
The Fulcrum, Vantage Way, Poole, Dorset BH12 4NU
tel: +44 (0) 1202 731555 fax: +44 (0) 1202 731666
info@ddcdolphin.co.uk
www.ddcdolphin.co.uk









## **DDC Pulpmatic Uno - Wall Recessed**

## Compact macerator for single medical pulp items with click & fit disinfection dosing system

The Pulpmatic Uno is a compact pulp macerator with capacity for 1-2 pulp items per cycle of 65 seconds, with the added benefit of our new click & fit disinfection dosing system. Developed for low demand areas of single bed wards and en-suite facilities, the Pulpmatic Uno is ideal for these situations, with its compact size, low noise and extremely low water and electricity consumption all housed in our new compact, sleek wall recessed model.

Click & Fit Dosing System - The new click & fit dosing system has been developed to make it even easier to ensure that sluice/dirty utility rooms are kept germ free. By simply opening the front flap and inserting the bottle you can leave the machine to accurately dose the EcoCleanse+ disinfectant into the machine, cleansing the drum after each cycle to meet EN1276 and EN13697 standards eliminating spores of up to 99.999%.

Unique Hands Free Operation - The Pulpmatic Uno uses our unique totally hands free technology offering completely hands free operation by means of a foot cup at the base of the machine and an optical sensor to close and start the machine. This makes it not only more convenient, but also more hygienic to use than other macerators and greatly reduces the potential for cross contamination.

Efficient and Economical - The Pulpmatic Uno offers greater efficiency, made possible

by the unique design features including the seven-blade maceration disc allowing faster and quieter operation. Not only does this help you cope better with workload peaks, but it translates into significant operating cost savings in terms of staff time and of water and electrical consumption per cycle.

Easy service programme - Service support is provided by our own specialist and highly skilled engineers and we offer a choice of fully comprehensive service contracts focused on individual customer needs.







## Key **Features**

Unique hands free open and close operation

Click & fit automatic disinfection dosing system

Dual antimicrobial action

Wall recessed

Efficient and economical

**Environmentally** friendly

Reduced water consumption

Easy to use

Easy to install

**Easy service** programme

development and improvement the described in this brochure

General	
Stainless Steel grade:	304
External Dimensions: w x h x d	460 x 1644 x 480 mm (flap closed)
Cutout Dimensions: w x h x d	452 x 1635 x 480 mm
Weight:	Dry 65kg
Capacity:	1-2 pieces medical pulp
Operation cycle time:	65 seconds
Decibel rating:	58.9 dba
Electrical	
Voltage: 110V	60Hz dedicated 15A breaker (UL version)
Voltage: 230V	50Hz dedicated 13A breaker
Voltage: 230V	60Hz dedicated 13A breaker
Voltage: 415V	N/A
Cable connection: 230V	3 core sheathed x 1.5mm², 2m long
Cable connection: 415V	N/A
Circuit Protection:	30mA RCD, GFI, IEC60364 compliant
Rated power input:	0.8kW
Electrical consumption per cycle:	0.01kW
Steam generator power:	N/A
Motor power:	0.75kW
Pump power:	0.125kW
Water	
Water connection:	3/4" BSP male fittings required
Water pressure:	Cold water 3/4" 0.5-10 bar
Water flow rate (Min):	5 Litres per minute at 1 bar
Water tank capacity:	15 Litres
Water consumption per cycle:	11 Litres
Waste connection:	N.B. 50mm (2")
Overflow connection:	22mm
Disinfection	
Recommended disinfection detergents:	EcoCleanse+
Disinfection tank volume: (internal)	1Litre
Disinfection dosing per cycle:	13 Milliliters